

Lab-Tips® Large Sampling Swab

LTP125S

Cleanroom laundered double layer 100% polyester knit swab

The LTP125S large knitted polyester swab is an excellent choice for general cleaning surface sampling. Its unique semi-flexible internal paddle provides greater surface coverage and grooves

on the handle make for easy snapping. Its cleanroom laundered knitted head provide ultra-low NVR and ion levels.



Key Attributes

- Cleanroom laundered 100% continuous filament 2-ply polyester knit head
- Long polypropylene handle with a large semi-flexible paddle head
- Thermally bonded head
- Grooves on handle for easy snapping
- Laundered and packaged in Class 10 {ISO Class 4} cleanroom

Benefits

- Critically low particles, fibers, ions and extractables
- Excellent chemical resistance
- High abrasion resistance while soft and non-abrasive
- No contaminating adhesives and good sorbency
- Autoclavable
- Free from Silicone, Amides and DOP

Applications

- Sampling surfaces for cleaning validation in many industries

Other sampling swab

- LTP125TOC—For TOC cleaning validation

Physical Characteristics

Head Material	Knitted polyester	Handle Thickness	2.8mm {0.110"}
Head Width	12.7mm {0.50"}	Handle Length	100mm {3.94"}
Head Thickness	5mm {0.196"}	Total Swab Length	125mm {4.92"}
Head Length	25mm {0.984"}	Head Bond	Thermal
Handle Material	Polypropylene	Handle Text	Berkshire
Handle Width	5.5mm {0.216"}		

www.berkshire.eu.com

Contact: Tel 44 (0) 870 757 2877
enquiries@berkshire.uk.com

Americas	Tel 1 413 528 2602	Fax 1 413 528 2614	info@berkshire.com
Europe	Tel 44 (0) 870 757 2877	Fax 44 (0) 870 757 2878	enquiries@berkshire.uk.com
SE Asia	Tel 65 6252 4313	Fax 65 6252 4312	enquiries@berkshire.com.sg
Japan	Tel 81 3 5827 2380	Fax 81 35827 2382	master@berkshire.co.jp

Technical Data:

ATTRIBUTE		LTP125S	TEST METHOD
Non-Volatile Residue	IPA mg/swab	0.075	IEST-RP.CC004.3, Sec 7.1.2
Ions (µg/swab)	Chloride	0.11	IEST-RP.CC004.3, Sec 7.2.2
	Sulfate	0.061	
	Nitrate	ND*	
	Phosphate	ND*	
	Fluoride	ND*	
	Potassium	0.010	
	Calcium	0.082	
	Sodium	0.032	
	Magnesium	0.011	

Notes:

- Data shown are typical values
- Due to differences in test methods applied and equipment utilized by different swab manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- * ND* = below the detection limit of the equipment.

Order Information:

Product	Item Number	Size	Swabs/pk	Pks/cs
Large Knitted Polyester Sampling	LTP125S.25	125mm	20	25


Other Berkshire products



Wipers



Mop Systems




Face Masks



Glove Liners



Documentation Systems



Swabs